


Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	BARD® MARQUEE® Disposable Core Biopsy Instruments and Kits
Device Model	MQ1410
Serial No.	1356886
Lot No.	0001356886
Manufacturer	BARD GmbH
Country of Origin	USA
Reference	https://www.bd.com/documents/alerts/PI_Bard-Marquee-disposable-core-biopsy-instruments-and-kits-recall_CL_EN.pdf
Device picture	
Reason of Recall	<p>NHRA initiates this FSN due to specific lots of BARD® MARQUEE® Disposable Core Biopsy Instruments and Kits have a potential fault in the Penetration Depth Switch. A Penetration Depth Switch detachment may result in:</p> <ul style="list-style-type: none"> • prolongation of the biopsy procedure as the user obtains a replacement device. • undershooting the lesion unexpectedly, compromising the sample quality and potentially leading to a misdiagnosis. • if the penetration is overshoot, the device may pierce adjacent tissue and/or organs causing injury.
Action should be taken	Please stop using the above mentioned defected medical device and contact the authorized representative Bahrain Pharmacy at bahpharm@batelco.com.bh to take necessary action for recall.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh