

## Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	CALiaGold® Calibrator CALiaGold® Control Set CALiaGold® Sample Diluent
Device Model	1151100 1151200 1151400
Lot No.	20137 20282 20360 20361 20362 20462 20853
Manufacturer	Sentinel.CH.SpA
Country of Origin	Italy
Reference	<a href="#">Link</a>
Reason of Recall	NHRA initiates this FSN due to an instability of the materials caused by some detergent lots used in the manufacturing of CALiaGold®
Action should be taken	Please stop using the above mentioned medical device and contact the authorized representative General Medical W.L.L at <a href="mailto:registration.medics@intercol.com">registration.medics@intercol.com</a> to take the necessary action for recall.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact [Medical\\_Devices@nhra.bh](mailto:Medical_Devices@nhra.bh)