

Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	ChemoLock Port
Device Model	CL-12 CL-70 CL2100
Lot No.	5819981 5842228 5816334
Manufacturer	ICU Medical
Country of Origin	USA
Reference	https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfres/res.cfm?id=195036
Reason of Recall	NHRA initiates this FSN due to a potential for unintended disconnection, inability to connect to or fully engage with injectors, due to a variation in the spring inside the Port, which may cause delay of therapy and exposure to caustic substances
Action should be taken	Please stop using the above mentioned medical device and contact the manufacturer to take the necessary action for recall.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh