

## Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details																											
<b>Device Name</b>	<b>Exmoor Plastics Nasal Splints</b>																										
<b>Lot No.</b>	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 60%;">Product</th> <th style="width: 40%;">Lots</th> </tr> </thead> <tbody> <tr> <td>Brown Nasal splint - REF N3</td> <td>70305 / 70653 / 71973</td> </tr> <tr> <td>Mackay Nasal splint - REF N4</td> <td>71237 / 71977</td> </tr> <tr> <td>Grimaldi Nasal splint(small) – REF N5</td> <td>70058 / 70662 / 71978</td> </tr> <tr> <td>Grimaldi Nasal splint(large) – REF N6</td> <td>70671 / 71871</td> </tr> <tr> <td>Shah Nasal splint (60mm x 27mm x 0.85mm) -REF N7</td> <td>70650 / 71996</td> </tr> <tr> <td>Shah Nasal splint (65mm x 31mm x 0.85mm) -REF N8</td> <td>71627 / 71238 / 71974 / 72320</td> </tr> <tr> <td>Shah Nasal splint (70mm x 37mm x 0.85mm) -REF N9</td> <td>71239 / 71873 / 70905</td> </tr> <tr> <td>Shah Nasal splint (60mm x 27mm x 1.25mm) -REF N10</td> <td>70300 / 72322</td> </tr> <tr> <td>Shah Nasal splint (65mm x 31mm x 1.25mm) -REF N11</td> <td>71626 / 71997</td> </tr> <tr> <td>Shah Nasal splint (70mm x 37mm x 1.25mm) -REF N12</td> <td>70151 / 71979</td> </tr> <tr> <td>Shah Nasal splint (60mm x 27mm x 1.75mm) -REF N13</td> <td>70367</td> </tr> <tr> <td>Nasal splint pairs – REF 23</td> <td>70306</td> </tr> </tbody> </table>	Product	Lots	Brown Nasal splint - REF N3	70305 / 70653 / 71973	Mackay Nasal splint - REF N4	71237 / 71977	Grimaldi Nasal splint(small) – REF N5	70058 / 70662 / 71978	Grimaldi Nasal splint(large) – REF N6	70671 / 71871	Shah Nasal splint (60mm x 27mm x 0.85mm) -REF N7	70650 / 71996	Shah Nasal splint (65mm x 31mm x 0.85mm) -REF N8	71627 / 71238 / 71974 / 72320	Shah Nasal splint (70mm x 37mm x 0.85mm) -REF N9	71239 / 71873 / 70905	Shah Nasal splint (60mm x 27mm x 1.25mm) -REF N10	70300 / 72322	Shah Nasal splint (65mm x 31mm x 1.25mm) -REF N11	71626 / 71997	Shah Nasal splint (70mm x 37mm x 1.25mm) -REF N12	70151 / 71979	Shah Nasal splint (60mm x 27mm x 1.75mm) -REF N13	70367	Nasal splint pairs – REF 23	70306
	Product	Lots																									
	Brown Nasal splint - REF N3	70305 / 70653 / 71973																									
	Mackay Nasal splint - REF N4	71237 / 71977																									
	Grimaldi Nasal splint(small) – REF N5	70058 / 70662 / 71978																									
	Grimaldi Nasal splint(large) – REF N6	70671 / 71871																									
	Shah Nasal splint (60mm x 27mm x 0.85mm) -REF N7	70650 / 71996																									
	Shah Nasal splint (65mm x 31mm x 0.85mm) -REF N8	71627 / 71238 / 71974 / 72320																									
	Shah Nasal splint (70mm x 37mm x 0.85mm) -REF N9	71239 / 71873 / 70905																									
	Shah Nasal splint (60mm x 27mm x 1.25mm) -REF N10	70300 / 72322																									
	Shah Nasal splint (65mm x 31mm x 1.25mm) -REF N11	71626 / 71997																									
	Shah Nasal splint (70mm x 37mm x 1.25mm) -REF N12	70151 / 71979																									
	Shah Nasal splint (60mm x 27mm x 1.75mm) -REF N13	70367																									
	Nasal splint pairs – REF 23	70306																									
<b>Manufacturer</b>	<b>Exmoor Plastics</b>																										
<b>Country of Origin</b>	UK																										
<b>Reference</b>	<a href="https://www.dha.gov.ae/190f288c-1cf7-4121-9b93-cf2f90ba605f">https://www.dha.gov.ae/190f288c-1cf7-4121-9b93-cf2f90ba605f</a>																										
<b>Reason of Recall</b>	NHRA initiates this FSN due to Due to a potential manufacturing defect concerning the nasal splints produced by the company, as Product coatings used in the production process have not been removed prior to packaging																										
<b>Action should be taken</b>	Please stop using the above mentioned medical device and contact the authorized representative to take the necessary action for recall.																										

**Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.**

For more information please contact [Medical\\_Devices@nhra.bh](mailto:Medical_Devices@nhra.bh)