


## Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	KN 95 Respirators
Device Model	All
Lot No.	All
Manufacturer	Jiangxi Benguo Technology Co.,Ltd.
Country of Origin	China
Reference	NHRA
Device picture	
Reason of Alert	NHRA initiates this FSN because the quality assurance certificate of the mentioned product is no longer recognized and therefore the quality, efficacy of the product and compliance to the international standards cannot be guaranteed. <b>(Not filtering 95% as per international standards).</b>
Action should be taken	The above mentioned mask can be used as a regular surgical mask and <b>Not</b> for COVID 19 frontline team. for more information, Please contact the distributor <b>Zad</b> at <a href="mailto:regulatory@groupzad.com">regulatory@groupzad.com</a> .

**Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.**

For more information please contact [Medical\\_Devices@nhra.bh](mailto:Medical_Devices@nhra.bh)