

Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the following medical device:

Device Details	
Device Name	Filter Tips, 1500 µl
Batch No.	0605020019, 0605020020
Manufacturer	Qiagen Sciences
Country of Origin	Germany
Reference	https://www.fda.gov/medical-devices/safety-communications/qiagen-recalls-filter-tips-use-qiasymphony-spas-instruments-due-potential-leak-which-may-result
Reason of Recall	NHRA initiates this FSN due to potential leak which may cause delayed or inaccurate results.
Action taken	In case of having the affected device, please stop using it and contact the local agent to take the necessary action.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh