


Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	Bard Max Core Disposable Core Biopsy Instrument
Lot No.	REDN0266
Manufacturer	Bard
Country of Origin	Switzerland
Reference	https://www.besse.com/BesseMedical/media/BesseMediaLibrary/pdfs/Recalls/BD-BIOPSY-PUNCH-RECALL.pdf?ext=.pdf
Device picture	
Reason of Recall	NHRA initiates this FSN due to the risk of having firing issues. An example of a firing issue is where the device does not completely fire, resulting in the inability to obtain tissue samples.
Action should be taken	In case of having the above defected medical device, please contact your authorized representative to take the necessary action for replacement.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh