


## Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	Footprint Ultra PK Suture Anchor
Device Model	72202902
Lot No.	2033391
Manufacturer	SMITH & NEPHEW
Country of Origin	COSTA RICA
Device picture	
Reason of Recall	NHRA initiates this FSN due to a packaging design related issue. The space between the device and the packaging permits excessive movement which could potentially lead to the displacement of the anchor's tip protector allowing the tip of the anchor to create a sterile breach.
Action should be taken	Please stop using the defected medical device and contact the authorized representative <b>Bahrain Pharmacy</b> at <a href="mailto:bahpharm@batelco.com.bh">bahpharm@batelco.com.bh</a> to take the necessary action for withdrawing.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact [Medical\\_Devices@nhra.bh](mailto:Medical_Devices@nhra.bh)