


## Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
<b>Device Name</b>	HEMASHIELD GOLD/HEMASHIELD PLATINUM
<b>Device Model</b>	M00202175828P0
<b>Lot No.</b>	18D04
<b>Serial No.</b>	1161644789 / 1161645334
<b>Manufacturer</b>	GETINGE
<b>Country of Origin</b>	France
<b>Reference</b>	<a href="https://www.healthycanadians.gc.ca/recall-alert-rappel-avis/hc-sc/2020/72811r-eng.php">https://www.healthycanadians.gc.ca/recall-alert-rappel-avis/hc-sc/2020/72811r-eng.php</a>
<b>Device picture</b>	
<b>Reason of Recall</b>	NHRA initiates this FSN due to intra-operative bleeding at the seam line of the branches of Hemashield grafts.
<b>Action should be taken</b>	Please stop using the above-mentioned defected medical device and contact the authorized representative to take the necessary action for recall.

**Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.**

For more information please contact [Medical\\_Devices@nhra.bh](mailto:Medical_Devices@nhra.bh)