


Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	R3 Acetabular Shells
Device Model	7133552
Lot No.	18HM16518
Manufacturer	SMITH & NEPHEW
Country of Origin	USA
Device picture	
Reason of Recall	NHRA initiated this FSN due to a manufacturing error. Use of the affected devices, which may not have been manufactured to their design specifications, may result in an intraoperative locking failure.
Action should be taken	Please stop using the above-mentioned defected medical device and contact the authorized representative Bahrain Pharmacy at pqc@bahrainpharmacy.com to take the necessary action for recall.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh