

Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	Neonatal Resuscitation System
Device Model	6430000
Lot No.	32054366 / 32055044 / 32056212 / 32013769
Manufacturer	Intersurgical LTD
Country of Origin	UK
Reference	http://www.salute.gov.it/imgs/C_17_AvvisiSicurezza_9943_azione_itemAzione0_files_itemFiles0_fileAzione.pdf
Reason of Recall	NHRA initiates this FSN due to incorrect or missing connection interfaces, this may prevent connection of the system to the particular equipment being used.
Action should be taken	Please stop using the above mentioned medical device and contact the authorized representative General Medicals at registration.medics@intercol.com to take the necessary action for replacement.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh