

Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	Seattle PAP plus Circuits, Coaxial Breathing Circuits, Pediatric Breathing Circuits, Pediatric Breathing Bag Set
Device Model	MP02260 MP07951 MP00351 MP00384
Manufacturer	Drager
Country of Origin	Germany
Reference	Link
Reason of Recall	NHRA initiates this FSN due to the adhesive joint did not conform to specification due to a manufacturing error. If system components become loose and/or detach from the hose during use, ventilation of the patient will be restricted.
Action should be taken	Please stop using the above mentioned medical device and contact the authorized representative Tariq Pharmacy at Registration.medics@intercol.com to take the necessary action for recall.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh