


Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	Pink Drop hand sanitizer gel 400ml
Lot No.	2020-0000154
Manufacturer	SAFIRES KOZMETIK
Country of Origin	Turkey
Reference	SFDA
Device picture	
Reason of Recall	The amount of sterile alcohol (Ethanol) is less than 60%, which is the minimum in hand sanitizers. Presence of methanol.
Action should be taken	Please stop using the above-mentioned defected product and contact the supplier to take the necessary action for recall

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh