

Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	DEFIGARD HD-7
Device Model	DEFIGARD HD-7
Manufacturer	Schiller Medical
Country of Origin	France
Reference	Link
Reason of Alert	NHRA initiates this FSN due to software faults in the DEFIGARD HD-7 monitor and defibrillator, which may cause incorrect timestamps or unexpected defibrillation modes, potentially delaying patient treatment.
Action should be taken	Please refer to "Actions to be taken by Customer/ User" in the attached FSN And for more information please contact the authorized representative General Medical W.L.L at registration.medics@intercol.com .

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh