

Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	Giraffe OmniBed Giraffe OmniBed Carestation
Affected Devices	All Giraffe OmniBed Carestations manufactured between January 2018 and November 2021 All Giraffe OmniBeds* and Giraffe OmniBed Carestations that have had heater doors replaced.
Manufacturer	GE Healthcare
Country of Origin	USA
Reference	Link
Reason of Alert	NHRA initiates this FSN due to improperly torqued screws securing heater doors in certain Giraffe OmniBed and Carestation devices, which may lead to loose doors, alarm activation, and potential patient injury if canopy movement is forced.
Action should be taken	Please refer to "Actions to be taken by Customer/ User" in the attached FSN And for more information please contact the authorized representative Gulf corporation for Technology W.L.L. at quality@gctbahrain.com .

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact Medical_Devices@nhra.bh